

HEMCO INDUSTRIES PTY LTD

OCCUPANT RESTRAINTS

ABN 12 004 746 856

7H3 Small Special Needs Harness For children from 14kg to 32kg

2 Feb 05

**This harness should only be used for children with special needs.
It is not a substitute for a general purpose child harness.**




Figure 1

General purpose child harnesses are not suitable for many children with special needs. The child is able to 'wriggle' out of the seat belt or, into a position such that the seat belt may become dangerous. The 7H3 is designed to envelop the passenger's trunk, thus following the wearer's movement and remaining in the correct, safest configuration. The wearer is unable to remove the harness. In exceptional cases, an extra 'anti-submarine' strap can be added to prevent the wearer sliding down & forward.

The harness is available for both cars and buses, in three sizes: small for children up to 32kg, and medium or large for children or adults over 32kg, with special needs.




Figure 2

Features

- It is often difficult to fit a harness onto a child with special needs. To make it easier for you, the 7H3 harness is designed so it can be fitted to the wearer before they get into the vehicle.
- Wide 48mm webbing has been used for increased comfort.
- Webbing in the neck area has been folded over so the edges do not chafe the wearer's neck.

Using the Harness

Fitting to the Wearer

Pass the open side of the harness over the wearer's right arm, with the shoulder straps above the arm and the chest strap below the arm. Pull the left shoulder strap over the wearer's head. Fasten the left side of the harness at the buckle as shown in Figure 3.




Figure 3

Adjustment

The harness is most effective when comfortably tight. All straps should be placed over bony areas, rather than soft tissue, the chest strap should sit across the centre of the breastbone, and the seat belt should sit over the wearer's hipbones. There are three sets of loops at the back to adjust the position of the chest strap. There are two sets of loops at the front to adjust the position of the seat belt. See Figures 2 & 5.

Attaching to vehicle

Place the child in position on the seat. Attach the snaphook on the tether strap to the anchor as shown in Figure 4. Pass the lap belt, or lap-sash combination (see over page) through the appropriate pair of loops at the front of the harness. Snap the tongue on the tether strap into the buckle on the harness. Adjust the seat belt so that it seats firmly across the wearer's hipbones, then tighten the harness using the tether strap adjustment, DO NOT over tighten such that the seat belt is raised off the child's hipbones.




Figure 4




Figure 5

HEMCO INDUSTRIES PTY LTD

OCCUPANT RESTRAINTS

ABN 12 004 746 856

7H3 Small Special Needs Harness For children from 14kg to 32kg

2 Feb 05

Converting Lap-sash Seat Belt to Lap

For lap-sash belts, fit the gated buckle to the belt so that the sash is doubled with the lap belt, as shown in Figure 6. Face the hollow side of the gated buckle towards the wearer. Position the gated buckle at the same height as the seat belt buckle on the other side of the wearer.




Figure 6

Installation of fittings

Vehicle with anchor point fitted

The Special Needs harness may be used with an existing child harness anchor if already fitted to the vehicle. Otherwise, remove the plastic keeper from the anchor point behind the seating position to be used. Remove the nut and reinforcing washer from the anchor bolt supplied. Screw the bolt and anchor into the anchor point thread. Do not over tighten (Max torque 20Nm).

Vehicle with no anchor point

For sedans only – Drill a 9mm hole in the parcel shelf on the centre line of the seating position as shown in Figure 7. Install the anchor bolt as shown in Figure 8. Do not over tighten (Max torque 20Nm).

For other vehicles – Refer to Hemco Industries at the above address for correct fitting to these vehicles, or for alterations needed for non-standard fitment.

WARNING: Restraint must not be attached to unsound metal or to wood, or synthetic structures. When drilling the 9mm hole do not drill through fuel lines, fuel tanks, electrical wiring or radio speakers.

Extras available from Hemco:

Anti-submarine strap, Fittings kit, Tether Extension Strap

WARNING:

•This accessory must be attached to a child restraint anchorage point in the car. Refer to vehicle Owner's Handbook for location of anchorage points. Refer to Hemco Industries at the above address if no anchorage points are fitted in the vehicle.

•Connect the restraint exactly as shown in the manufacturer's instructions.

•Should any of these parts become worn, frayed, contaminated, damaged, or if the restraint has been in a severe crash, the device, including the fitting kits, should be completely destroyed, even if no damage is obvious.

•Tether and extension straps may be damaged by unsecured cargo in a collision. Hazardous cargo should be secured or removed before the restraint is used.

•Ensure the webbing does not pass over sharp edges.

•Clean with mild soap and water. Avoid contamination of webbing with polishes, oils and other chemicals.

• Do not alter harness (consult Hemco Industries at the above address)




Figure 7




Figure 8

List of items included in this kit:

- 1 7H3 Small Harness
- 1 Adjustable Tether Strap
- 1 Gated Buckle
- 1 Hex Head Set Screw, 30mm x 5/16" - 18 UNC
- 1 Spring Washer 5/16"
- 1 Spacer, 5mm
- 1 Spacer, 10mm
- 1 Reinforcing Plate
- 1 Nut, 5/16" - 18 UNC

Contact Hemco Industries at the above address for replacement items or optional extras.